



Recombinant *Schizosaccharomyces pombe* SWI/SNF and RSC complexes subunit *ssr3* (*ssr3*)

Product Code	CSB-MP879135SXV
Abbreviation	<i>ssr3</i>
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P7S3
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSNNSRLPEN GVQSGNGEDA ELKKS MRIIE REIPDSLLEK IPEAEDYIAL QDLERKLDL IVRKRFDLQD SLSRNSHKTR ILRMYIHSTV ANQSWQQKGE NQENNSGDIN SLPIPEWTLH IEGRLLVNP DEDDKAFELA PFTNFFRKIA IQILRSDDLY PSGNYVEW NK LPDNSNTSNG ITVTRKGDQS VDVKIMLYPE EHPERYKLSK AFANILGIRE GTRPDIVSYL WQYIKFHRLQ DMEEKRLINC DKALRDLFEA DRLYFPRIPE LMNRFLEPID PIVIPYTINV SEHTVEKVTI FDIRINTEDP RHSQIRSFLA TMMSQDKIRS IDDKLTEL IQ AITYSQSKYD FMKKFSESPI EFINEWIESQ SRDLEIVLDG TNMNYAEKRS ADYYQQPWVH ESAFHYLNLL NSKKQQSVLN ASAKK
Source	Mammalian cell
Target Names	<i>ssr3</i>
Protein Names	Recommended name: SWI/SNF and RSC complexes subunit <i>ssr3</i>
Expression Region	1-425
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.