



Recombinant Transcription termination factor Rho (rho)

Product Code	CSB-BP879155XBL
Abbreviation	rho
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9PA21
Product Type	Recombinant Protein
Immunogen Species	Xylella fastidiosa (strain 9a5c)
Purity	>85% (SDS-PAGE)
Sequence	MPAHKLEIA EQLNIHEGVA RARKQDVIFA LLKVLTRHGE GVFADGVLEI LPDGFGLRA AEASYLAGPD DTYISPSQIR RFNLRTGDHL SGRIRFPKDG ERYFALSIVD TINGEPLEAS KNKILFENLT PLFPRKRFRL ERADGSTEDI TGRILDLMAP QGKGQRALIV SPPKAGKTM LMQQVATAIT NHPEVHLIVL LIDERPEEVT EMQRTVRGEV ISSTFDEPAA RHVQVAEMVI ERAKRLIEHK KDVVILLDSI TRLARAYNNV VPSSGKVL TGVDANALHRP KRFFGAARNV EEGSLTIIA TALIDTGSKM DEVIYEEFKG TGNCELHLNR RIAEKRVYPA IDINRSGTRR EDLLIESDLL QKIWILRKLL HPMDEIAAIE FVLDMKMKNTK SNDEFFGSMK R
Source	Baculovirus
Target Names	rho
Protein Names	Recommended name: Transcription termination factor Rho EC= 3.6.4.- Alternative name(s): ATP-dependent helicase Rho
Expression Region	1-411
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.