



Recombinant Rat Cyclin-L1 (Ccn1)

Product Code	CSB-MP879508RA
Abbreviation	Ccn1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9R1Q2
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MASGPHPTST AAAASASSAA PSAGGSSSGT TTTTTTTTGG ILIGDRLYSE VSLTIDHSVI PEERLSPTPS MQDGLDLPSE TDLRILGCEL IQAAGILLRL PQVAMATGQV LFHRFFYSKS FVKHSFEIVA MACINLASKI EEAPRRIRDV INVFHHLRQL R GKRTPSPLI LDQNYINTKN QVIKAERRVL KELGFCVHVK HPHKIIVMYL QVLECERNQT LVQTAWNYMN DSLRTNVFVR FQPETIACAC IYLAARALQI PLPTRPHWFL LFGTTEEEIQ EICIETLRLY TRKKPNYELL EKEVEKRKVA LQEAKLKAKG LNLDTGTPALS TLGGFSPASK PSSPREVKAE EKSPVSINVK TVKKEPEDRQ QASKSPYNGV RKDSKRSRNS RSASRSRSRT RSRSRSHTPR RHYNNRRSRS GTYSSRSRSR SRSHSESPRR HHNHGSPHLK AKHTREDLKS SNRHGHKRKK SRSRSQSKTR DHSVDTKKHR HERGHHRDRR ERSRSFERSH KGKHHGGSRS GHGRHRR
Source	Mammalian cell
Target Names	Ccn1
Protein Names	Recommended name: Cyclin-L1 Short name= Cyclin-L Alternative name(s): Cyclin Ania-6a
Expression Region	1-527
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.