



Recombinant Escherichia coli Uncharacterized protein yubM (yubM)

Product Code	CSB-MP879667ENV
Abbreviation	yubM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9S4W2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSVTESKAKT ERKSSRPAK TQETVLSALL AQTEEVSVPL ASLIKSPLNV RTVPYSAESV SELADSIKGV GLLQNLVVHA LPGDRYGVAA GRRRLAALNM LAERDIIQVD WPVRVKVIPQ ELATAASMTE NGHRRDMHPA EQIAGFRAMA QEGKTPAQIG DLLGYSPRHV QRMLKLADLA PVILDALAED RITTEHCQAL ALENDTARQV QVFEEACQSG WGGKPDVRVI RNLITSEVA VKDNTKFRFV GADAFSPDEL RTDLFSDDEG GYVDCVALDA ALLEKLQAVA EHLREAEGWE WCAGRMEPVG ECREDSRAYR NLPEPEAVLT EEEEEERLNEL MMRYDALENQ CEESDLLAAE MKLIDCMAKV RAWTPEMRAG SGVVVSWRYG NVCVQRGVQL RSEDDVTDDA DRTEQVQEKA SVEEISLPLL TKMSSERTLA VQAALMQQPD KSLALLTWTL CLNVFGSGAY SKPAQISLEC KHYSLTSDAP SGKEGA AFMA LMAEKARLVV LLPEGWSRDM TTFLSLSQEV LLSLLSFCTA CSIHGVQTRE CGHTSR SPLD PLETAIGFHM RDWWQPTKAN FFGHLKPKQI IAALNEAGLS GAARDAEKMK KGDAAEHAEF HMKDNRWVPG WMCAPRPQTD ATERTDNLAD AA
Source	Mammalian cell
Target Names	yubM
Protein Names	Recommended name: Uncharacterized protein yubM
Expression Region	1-652
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.