



Recombinant Arabidopsis thaliana Putative F-box/LRR-repeat protein 23 (FBL23)

Product Code	CSB-YP879689DOA
Abbreviation	FBL23
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9S9V9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MAPHRSRQRN PRFKLWRYDD EPLAIAETMP ELRHLKLFNG GLTNLRLEAI LDNCVHLVHL DLRRCFNINL LGDLEKRCSE RIRDLRRPDD STADSPFDAS SDIYSAGEDD YDFYSDDSDG TLDHFEECK TVSVATGNES VYVLSFVRS NFIYPKKIAS SSSLSSTPLP SLMKDKKLRN WAELPSKLT SILLRLGAIE ILQNAQKVCK PWHRVCKDPS MWRKIDIDNR NDRAAFKYDL ESMCRHAVDR SHGGLIEIEI WYYGTNDLIM YIADRSSNLK SLGLVRCFPI TDEGVAKAVS KVPLLEYLEV SYCLFSGESL RDIGRSCP NLKTLKLNRAPE IMFSNSGFDD NAKAIAESMP ELRHLQLLGN GLTNKGLNAI LDGCPHLEHL DLRQCFNINL VGDLKKRCFE RIKDLRCPND SDDSDDGSD NGAVTTFGSQ FNETDYHWY
Source	Yeast
Target Names	FBL23
Protein Names	Recommended name: Putative F-box/LRR-repeat protein 23
Expression Region	1-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.