



# Recombinant Arabidopsis thaliana Probable protein phosphatase 2C 6 (At1g16220)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-YP879691DOA  |
| <b>Abbreviation</b>      | At1g16220  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q9SA22   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Arabidopsis thaliana (Mouse-ear cress)   |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MGLCHSKIDK TTRKETGATS TATTTVERQS SGRLRRPRDL YSGGEISEIQ<br>QVVGRLVGNG SSEIACLYTQ QGKKG TNQDA MLVWENFCSR SDTVLCGVFD<br>GHGPF GHMVS KRVRDMLPFT LSTQLKTTSG TEQSSSKNGL NSAPTCVDEE<br>QWCELQLCEK DEKLFPEMYL PLKRALLKTC QQMDKELKMH PTINCFCSGT<br>TSVTVIKQGK DLVVG NIGDS RAVLATRDQD NALVAVQLTI DLKPDLPSES<br>ARIHRCKGRV FALQDEPEVA RVWLPNSDSP GLAMARAFGD FCLKDYGLIS<br>VPDINYHRLT ERDQYIILAT DGVWDVLSNK EAVDIVASAP SRDTAARAVV<br>DTAVRAWRLK YPTSKNDDCA VVCLFLEDTS AGGTVEVSET VNHSHEESTE<br>SVTITSSKDA DKKEEASTET NETVPVWEIK EEKTPESCRI ESKKTTLAEC<br>ISVKDDEEWS ALEGLTRVNS LLSIPRFFSG ELRSSSWRKW L |
| <b>Source</b>            | Yeast  |
| <b>Target Names</b>      | At1g16220  |
| <b>Protein Names</b>     | Recommended name: Probable protein phosphatase 2C 6 Short name=<br>AtPP2C06 EC= 3.1.3.16   |
| <b>Expression Region</b> | 1-491  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life  |



of lyophilized form is 12 months at -20°C/-80°C.