



Recombinant Arabidopsis thaliana Zinc finger CCCH domain-containing protein 34 (At3g06410)

Product Code	CSB-MP879807DOA
Abbreviation	At3g06410
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SQU4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MERYGRPGEE GSRSDPSLEW TSHGGETAVE APMWRLGLSG GGGGGGGGES YPERPDEPDC IYYLRTGVCG YGSRCRFNHP RDRGAVIGGV RGEAGALPER MGHPVCQHFM RTGTCKFGAS CKYHHPRQGG GGGSVAPVSL SYLGYPLRPG EKECSYYLRT GQCKFGLTCR FNHPVPLAVQ GPPQQPQQQQ PQPQPQLQTI YPTLQSQSIP SSQQYGLVLT RPSFLTGSYL QSPYGPPMVL PPGMVPYSGW NPYQASLSAM PSPGTQPSIG SSSIYGLTPL SPSATAYTGT YQSVSSNST SKEFPQRPDQ PECQYFMRTG DCKFGSSCRY HHPVDAVPPK TGIVLSSIGL PLRPGVAQCT HFAQHGICKF GPACKFDHSM SSSLSYSPSA SSLTDMPVAP YPIGSSSLSG SSAPVSSSNE PTKEAVTPAV SSMVSGLSRP EPAETSGDSA SVSGSIEAKT SS
Source	Mammalian cell
Target Names	At3g06410
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 34 Short name= AtC3H34 Alternative name(s): Zinc finger CCCH domain-containing protein ZFN-like 2
Expression Region	1-462
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.