



Recombinant Arabidopsis thaliana Putative F-box protein At3g10240 (At3g10240)

Product Code	CSB-MP879822DOA
Abbreviation	At3g10240
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SS35
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEQQEEKRKI KAYQRKSKRS KSGSSSIPLD LVSEILLRLP EKSVARFRCV SKPWSSITTE PYFINLLTTR SPRLLLCFKA NEKFFVSSIP QHRQTFETWN KSHSYSQLID RYHMEFSEEM NYFPPTESVN GLICFQESAR LIVWNPSTRQ LLILPKPNGN SNDLTIFLG Y DPVEGKHKVM CMEFSATYDT CRVLTGSAQ KLWRTVKTHN KHRSDYYDSG RCINGVVYHI AYVKDMCVWV LMSFDVRSEI FDMIPLSSD VH KDVLIDYN GRLACVGREI IEKNGIRLWI LEKHNKWSSK DFLAPLVHID KSLSTNKFL KGFTHAGEII YVESMFHKSA KIFFYDPVRN TSRRFELKGF TDDEFVLSNE HGYTLHVFPN HVESQISFT
Source	Mammalian cell
Target Names	At3g10240
Protein Names	Recommended name: Putative F-box protein At3g10240
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.