



# Recombinant Arabidopsis thaliana High mobility group B protein 6 (HMGB6)

<b>Product Code</b>	CSB-BP879848DOA
<b>Abbreviation</b>	HMGB6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9SUP7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MATNADPAPT KKPRNSRKAL KQKNELVETP PSPVSVKGKS AKSFEQDLME MQTMLEKMKI EKDKTEELLK EKDEILRKKE EELETRDAEQ EKLKVELKKL QKMKEFKPNM TFACGQSSLT QAEQEKANKK KKKDCPETKR PSSSYVLWCK DQWTEVKKEN PEADFKETSN ILGAKWKSLA AEDKKPYEER YQVEKEAYLQ VIAKEKREKE AMKLEDDQK QRTAMELLDQ YLNFVQAEQ DNKKKNKKEK DPLKPKHPVS AFLVYANERR AALREENKSV VEVAKITGEE WKNLSDKKA PYEKVAKKNK ETYLQAMEEY KRTKEEEALS QKKEEEELLK LHKQEALQML KKKEKTDNLI KKEKATKAAA NENVDPNPKP KPASSYFLFS KDERKKLTEE RPGTNNATVT ALISLKWKEL SSEEKQVYNG KAAKLMEAYK KEVEAYNKKS AATTSS
<b>Source</b>	Baculovirus
<b>Target Names</b>	HMGB6
<b>Protein Names</b>	Recommended name: High mobility group B protein 6 Alternative name(s): Nucleosome/chromatin assembly factor group D 06 Short name= Nucleosome/chromatin assembly factor group D 6 WRKY transcription factor 53 Short name= AtWRKY53
<b>Expression Region</b>	1-456
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.