



Recombinant Arabidopsis thaliana E3 ubiquitin-protein ligase RHA1B (RHA1B)

Product Code	CSB-EP879850DOA-B
Abbreviation	RHA1B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SUS5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGLPTDFKEL QIPGYVLKTL YVIGFFRDMV DALCPYIGLP SFLDHNETS SDPTLALST SATLANELIP VVRFSDLLTD PEDCCTVCLS DFVSDDKIRQ LPKCGHVFHH RCLDRWIVDC NKITCPICRN RFLPEEKSTP FDWGTSDWFR DEVESTN
Source	E.coli
Target Names	RHA1B
Protein Names	Recommended name: E3 ubiquitin-protein ligase RHA1B EC= 6.3.2.- Alternative name(s): RING-H2 zinc finger protein RHA1b
Expression Region	1-157
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.