



Recombinant Arabidopsis thaliana Putative pentatricopeptide repeat-containing protein At1g10330 (PCMP-E71)

| | |
|--------------------------|---|
| Product Code | CSB-BP879879DOA |
| Abbreviation | PCMP-E71 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9SY75 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MRCTSFSSFS LSLEDALHLL QRFLYSSNQI KQIHTVLLTS NALVASRWKT KCVYNTLIRS YLTTGEYKTS LALFTHMLAS HVQPNNLTFP SLIKAACSSF SVSYGVALHG QALKRGFLWD PFVQTSFVRF YGEVGDLESS RKMFDLILNP CVVACNSLLD ACGRNGEMDY AFEYFQRMPV TDVVSWTTVI NGFSKKGLHA KALMVFGEMI QNERAVITPN EATFVSVLSS CANFDQGGIR LGKQIHGYVM SKEIILTTTL GTALLDMYGK AGDLEMALTI FDQIRDKKVC Awnaiisala SNGRPKQALE MFEMMKSSYV HPNGITLLAI LTACARSKLV DLGIQLFSSI CSEYKIIPTS EHYGCVVDLI GRAGLLVDAA NFIQSLPFEP DASVLGALLG ACKIHENTEL GNTVGKQLIG LQPQHCGQYV ALSTFNALDS NWSEAEKMRK AMIEAGIRKI PAYSVLT |
| Source | Baculovirus |
| Target Names | PCMP-E71 |
| Protein Names | Recommended name: Putative pentatricopeptide repeat-containing protein At1g10330 |
| Expression Region | 1-467 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.