



Recombinant Arabidopsis thaliana BTB/POZ and TAZ domain-containing protein 3 (BT3)

Product Code	CSB-BP879883DOA
Abbreviation	BT3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SYL0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSSSTKNIPK PPPLPCITYQ RFQSSTRKPS SLMRLVPKEA LETWDKLFKE GSGADTYVET DNKSHFPAHS SVLAAASPVI ATLLNQSRDK NGNTYLKIHG VPCEAVYMFY RFLYSSCYEE EEMKKFVLHL LVLSHCYSPV SLKRLCVEL DQGWINKENV IDVLQLARNC DVTRICFVCL SMVIKDFKSV SSTEHWKVMK RSNPLEQEL IEAVIESDSR KQERRRLEE REVYLQLYEA MEALVHICRE GCGTIGPRDK ALKGSHTVCK FPACKGLEGA LRHFLGCKSR ASCSHCKRMW QLLQLHSCIC DDSNSCKVSL CWNFKEKMKK LSKKEQSKWR LLVENIIRAR NSLGPFFSSRS SGLI
Source	Baculovirus
Target Names	BT3
Protein Names	Recommended name: BTB/POZ and TAZ domain-containing protein 3 Alternative name(s): BTB and TAZ domain protein 3
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.