



# Recombinant Arabidopsis thaliana GDSL esterase/lipase At4g30140 (At4g30140)

<b>Product Code</b>	CSB-EP879895DOA-B
<b>Abbreviation</b>	At4g30140
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9SZW7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	QQ TPCYFVFGDS VFDNGNNAL NTKAKVNYLP YGIDYFQGPT GRFSNGRNIP DVIAELAGFN NPIPPFAGAS QAQANIGLNY ASGAGGIREE TSENMGERIS LRQQVNNHFS AITAAVPLS RLRQCLYTIN IGSNDYLNYY FLSPPTLARR LFNPDQYARS LISLYRIYLT QLYVLGARNV ALFGIGKIGC TPRIVATLGG GTGCAEEVNQ AVIIFNTKLLK ALVTDFNPKP GAMFTYVDLF SGNAEDFAAL GITVGDRSCC TVNPGHEELCA ANGPVCPDRN KFIFWDNVHT TEVINTVVAN AAFNGPIASP FNISQLVN
<b>Source</b>	E.coli
<b>Target Names</b>	At4g30140
<b>Protein Names</b>	Recommended name: GDSL esterase/lipase At4g30140 EC= 3.1.1.- Alternative name(s): Extracellular lipase At4g30140
<b>Expression Region</b>	29-348
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.