



Recombinant Acyrthosiphon pisum secondary endosymbiont phage 1 Putative protein p9 (9)

Product Code	CSB-EP879912AKAV-B
Abbreviation	9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9T1T9
Product Type	Recombinant Protein
Immunogen Species	Acyrthosiphon pisum secondary endosymbiont phage 1 (Bacteriophage APSE-1)
Purity	>85% (SDS-PAGE)
Sequence	MLNIKCIPRK ALVLVLGGGM ILPNISPAMI NDTAGCAGAG ISYPYIAGAN DNPLFISSDH RPTNNINQVM MMEIIRTINN QTTFQTVRAA RHHIIPYQVL RDFYNAVITR GQSNPLELRT FGQILLHLVD NASRLYQNTI SLSQLDDARS VALDFMHGVV LGDETRFSEQ QVSGLYAIR LFTWF PANIF IGPDTSYRFD DPGNDFEUNA RSIVGAENFN RLEALHNL MV AYINKRNTLD HGWSVIMLAR MVSMVSLFKA LDNVVTISNA QNWFEPDPWY GDERVGIGFN GRVEYVDVQN LNPPENMRKK RENKNNFIVK RLPTIADPYC LMAPKIVKKK LVASLFTDIN PDLNTNKEPK HDEL RK
Source	E.coli
Target Names	9
Protein Names	Recommended name: Putative protein p9
Expression Region	1-366
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.