



Recombinant Human Phosphatidylserine decarboxylase proenzyme, mitochondrial (PISD)

Product Code	CSB-EP880077HU-B
Abbreviation	PISD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UG56
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MMQCSEARQGPELRAAKWLHFPQLALRRRLGQLSCMSRPALKLRSWPLTVL YYLLPFGALRPLSRVGRPVSRVALYKSVPTRLLSRAWGRLNQVELPHWLRR PVYSLYIWTFGVNMKEAAVEDLHHRNLSEFFRRKLKQPQARVPCGLHSHVISPS DGRILNFGQVKNCEVEQVKGVTYSLESFLGPRMCTEDLPFPPAASCDSFKNQ LVTREGNELYHCVIYLAPGDYHCFHSPDWTVSHRRHFPGLMSVNPGMAR WIKELFCHNERVVLTDGDWKHGFFSLTAVGATNVGSIRIYFDRDLHTNSPRHSK GSYNDFSFTHTNREGVPMRKGEHLGEFNLGSTIVLIFEAPKDFNFQLKTGQKI RFGEALGSL
Source	E.coli
Target Names	PISD
Protein Names	Recommended name: Phosphatidylserine decarboxylase proenzyme EC= 4.1.1.65Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Isoform 2
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.