



Recombinant Human Nuclear receptor-interacting protein 2 (NRIP2)

Product Code	CSB-BP880083HU
Abbreviation	NRIP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BQI9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLFIFPLSLP WRPSCWKESC STGQRQAGRS REDSVTPPPS SPWPTPPAGA MSTKQEARRD EGEARTRGQE AQLRDRAHLS QQRRLKQATQ FLHKDSADLL PLDSLKRLGT SKDLQPRSVI QRRLVEGNPN WLQGEPPRMQ DLIHGQESRR KTSRTEIPAL LVNCKCQDQL LRVAVDTGTQ YNRISAGCLS RLGLEKRVLK ASAGDLAPGP PTQVEQLELQ LGQETVCSA QVVDAESPEF CLGLQTLLSL KCCIDLEHGV LRLKAPFSEL PFLPLYQEPG Q
Source	Baculovirus
Target Names	NRIP2
Protein Names	Recommended name: Nuclear receptor-interacting protein 2
Expression Region	1-281
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.