



# Recombinant Human Coronin-1B (CORO1B)

<b>Product Code</b>	CSB-BP880087HU
<b>Abbreviation</b>	CORO1B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BR76
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSFRKVVRRQS KFRHVFVGGQPV KNDQCYEDIR VSRVTWDSTF CAVNPKFLAV IVEASGGGAF LVLPLSKTGR IDKAYPTVCG HTGPVLDIDW CPHNDEVIAS GSEDCTVMVW QIPENGLTSP LTEPVVVLEG HTKRVGIIAW HPTARNVLLS AGCDNVVLIW NVGTAEELYR LDSLHPDLIY NVSWNHNGSL FCSACKDKSV RIIDPRRGTL VAEREKAHEG ARPMRAIFLA DGKVFTTGFS RMSERQLALW DPENLEEPMA LQELDSSNGA LLPFYDPDTS VVYVCGKGDGDS SIRYFEITEE PPYIHFLNTF TSKEPQRGMG SMPKRGLEVS KCEIARFYKL HERKCEPIVM TVPRKSDLFQ DDLYPDTAGP EAALEAEWV SGRDADPILI SLREAYVPSK QRDLKISRRN VLSDSRPAMA PGSSHLGAPA STTTAADATP SGSLARAGEA GKLEEVMQEL RALRALVKEQ GDRICRLEEQ LGRMENGDA
<b>Source</b>	Baculovirus
<b>Target Names</b>	CORO1B
<b>Protein Names</b>	Recommended name: Coronin-1B Alternative name(s): Coronin-2
<b>Expression Region</b>	1-489
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.