



Recombinant Human Coronin-1B (CORO1B)

Product Code	CSB-YP880087HU
Abbreviation	CORO1B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BR76
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSFRKVVQRQS KFRHVFVGGPV KNDQCYEDIR VSRVTWDSTF CAVNPKFLAV IVEASGGGAF LVLPLSKTGR IDKAYPTVCG HTGPVLDIDW CPHNDEVIAS GSEDCTVMVW QIPENGLTSP LTEPVVVLEG HTKRVGIIAW HPTARNVLLS AGCDNVVLIW NVGTAEELYS LDSLHPDLIY NVSWNHNGSL FCSACKDKSV RIIDPRRGTL VAEREKAHEG ARPMRAIFLA DGKVFTTGFS RMSERQLALW DPENLEEPMA LQELDSSNGA LLPFYDPDTS VVYVCGKGDGDS SIRYFEITEE PPYIHFLNTF TSKEPQRGMG SMPKRGLEVS KCEIARFYKL HERKCEPIVM TVPRKSDLFQ DDLYPDTAGP EAALEAEWV SGRDADPILI SLREAYVPSK QRDLKISRRN VLSDSRPAMA PGSSHLGAPA STTTAADATP SGSLARAGEA GKLEEVMQEL RALRALVKEQ GDRICRLEEQ LGRMENGDA
Source	Yeast
Target Names	CORO1B
Protein Names	Recommended name: Coronin-1B Alternative name(s): Coronin-2
Expression Region	1-489
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.