



Recombinant Human Protein NipSnap homolog 3B (NIPSNAP3B)

Product Code	CSB-MP880093HU
Abbreviation	NIPSNAP3B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BS92
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLVLRSGGTK ALASRTLAPQ VCSSFATGPR QYDGTIFYEFR TYYLKPSNMN AFMENLKKNI HLRTSYSELV GFWSVEFGGR TNKVFHIWKY DNFAHRAEVR KALANCKEWQ EQSIIPNLAR IDKQETEITY LIPWSKLEKP PKEGVYELAV FQMKPGGPAL WGDAFERAIN AHVNLGYTKV VGVFHTTEYGE LNRVHVLWWN ESADSRAAGR HKSHEDPRVV AAVRESVNYL VSQQNMLLIP ASFSPLK
Source	Mammalian cell
Target Names	NIPSNAP3B
Protein Names	Recommended name: Protein NipSnap homolog 3B Short name= NipSnap3B Alternative name(s): SNAP1
Expression Region	1-247
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.