



Recombinant Human Protease-associated domain-containing protein 1 (PRADC1)

Product Code	CSB-MP880094HU
Abbreviation	PRADC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BSG0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	HGFRIHDYL YFQVLSPGDI RYIFTATPAK DFGGIFHTRY EQIHLVPAEP PEACGELSNG FFIQDQIALV ERGGCSFLSK TRVVQEHGGR AVIISDNAV NDSFYVEMIQ DSTQRTADIP ALFLLGRDGY MIRRSLEQHG LPWAIISIPV NVTSIPTFEL LQPPWTFW
Source	Mammalian cell
Target Names	PRADC1
Protein Names	Recommended name: Protease-associated domain-containing protein 1 Alternative name(s): Protease-associated domain-containing protein of 21 kDa Short name= hPAP21
Expression Region	22-188
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.