



Recombinant Human Mini-chromosome maintenance complex-binding protein (MCMBP)

Product Code	CSB-MP880100HU
Abbreviation	MCMBP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BTE3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPCGEDWLSH PLGIVQGFFA QNGVNPDWEK KVIEYFKEKL KENNAPKWVP SLNEVPLHYL KPNSFVKFRC MIQDMFDPEF YMGVYETVNQ NTKAHLHFG KYRDVAECGP QQELDLNSPR NTTLERQTFY CVPVPGESTW VKEAYVNANQ ARVSPSTSYT PSRHKRSYED DDDMDLQPNK QKDQHAGARQ AGSVGGLQWC GEPKRLETEA STGQQNLNLSN LSSPFDLNFN LPGEKGPACL VKVYEDWDCF KVNIDILELYG ILSVDPVLSI LNNDERDASA LLDPMECTDT AEEQRVHSPP ASLVPRIHVI LAQKLQHINP LLPACLNKEE SKTCKFVSSF MSELSPVRAE LLGFLTHALL GDSLAAEYLI LHLISTVYTR RDVLPLGKFT VNLSGCPRNS TFTEHLYRII QHLVPASFRL QMTIENMNHL KFIPHKDYTA NRLVSGLLQL PSNTSLVIDE TLLEQGQLDT PGVHNVTALS NLITWQKVDY DFSYHQMEFP CNINVFITSE GRSLLPADCQ IHLQPQLIPP NMEEYMNSSL SAVLPSVLNK FRIYLTLLRF LEYSISDEIT KAVEDDFVEM RKNDPQSITA DDLHQLLVVA RCLSLSAGQT TLSRERWLRA KQLESLRRTR LQQQKCVNGN EL
Source	Mammalian cell
Target Names	MCMBP
Protein Names	Recommended name: Mini-chromosome maintenance complex-binding protein Short name= MCM-BP Short name= MCM-binding protein
Expression Region	1-642
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.