



Recombinant Human Uncharacterized protein C1orf50 (C1orf50)

Product Code	CSB-BP880109HU
Abbreviation	C1orf50
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BV19
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEDAAAPGRT EGVLERQGAP PAAGQGGALV ELTPTPGGLA LVSPYHTHRA GDPLDLVALA EQVQKADEFI RANATNKLV IAEQIQHLQE QARKVLEDAH RDANLHHVAC NIVKKPGNIY YLYKRESGQQ YFSIISPKEW GTSCPHDFLG AYKLQHDLSW TPYEDIEKQD AKISMMDTLL SQSVALPPCT EPNFQGLTH
Source	Baculovirus
Target Names	C1orf50
Protein Names	Recommended name: Uncharacterized protein C1orf50
Expression Region	1-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.