



Recombinant Human Alpha N-terminal protein methyltransferase 1A (METTL11A)

Product Code	CSB-BP880111HU
Abbreviation	METTL11A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BV86
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTSEVIEDEK QFYSKAKTYW KQIPPTVDGM LGGYGHISSI DINSSRKFLQ RFLREGPNKT GTSCALDCGA GIGRITKRL LPLFREVD MV DITEDFLVQA KTYLGEEGKR VRNYFCCGLQ DFTPEPDSYD VIWIQWVIGH LTDQH LAEFL RRCKGSLRPN GIIVIKDNMA QEGVILDDVD SSVCRDL DVV RRIICSAGLS LLAEERQENL PDEIYHVYSF ALR
Source	Baculovirus
Target Names	NTMT1
Protein Names	Recommended name: Alpha N-terminal protein methyltransferase 1A EC=2.1.1.n5 Alternative name(s): Methyltransferase-like protein 11A N-terminal RCC1 methyltransferase X-Pro-Lys N-terminal protein methyltransferase 1A Short name=
Expression Region	1-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.