



Recombinant Human RNA-binding protein 4 (RBM4)

Product Code	CSB-YP880118HU
Abbreviation	RBM4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BWF3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MVKLFIGNLP REATEQEIRS LFEQYGKVL E CDIKNYGFV HIEDKTAAED AIRNLHHYKL HGVNINVEAS KNKSKTSTKL HVGNISPTCT NKELRAKFEE YGPVIECDIV KDYAFVHMER AEDAVEAIRG LDNTEFQGKR MHVQLSTSRL RTAPGMGDQS GCYRCGKEGH WSKECPIDRS GRVADLTEQY NEQYGAVRTP YTMSYGDSLY YNNAYGALDA YYKRCRAARS YEAVAAAAAS VYNYAEQTLS QLPQVQNTAM ASHLTSTSLD PYDRHLLPTS GAAATAAAAA AAAAAVTAAS TSYYGRDRSP LRRATAPVPT VEGEGYGYGHE SELSQASAAA RNSLYDMARY EREQYADRAR YSAF
Source	Yeast
Target Names	RBM4
Protein Names	Recommended name: RNA-binding protein 4 Alternative name(s): Lark homolog Short name= hLark RNA-binding motif protein 4 RNA-binding motif protein 4a
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.