



Recombinant Human Single-stranded DNA-binding protein 3 (SSBP3)

Product Code	CSB-BP880121HU
Abbreviation	SSBP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BWW4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFAKKGKGS AV PSDGQAREKL ALYVVEYLLH VQAQKSAQTF LSEIRWEKNI TLGEPGFLH SWWCVFWDLY CAAPERDTC EHSSEAKAFH DYSAAAAPSP VLGNIPPNDG MPGGPIPPGF FQGPPGSQPS PHAQPPPHNP SSMMGPHSQP FMSPRYAGGP RPIRMGNQP PGGVPGTQPL LPNSMDPTRQ QGHPNMGGS M QRMNPPRGMG PMGPGPQNYG SGMRPPPNL GPAMPGINMG PGAGRPWPNP NSANSIPYSS SSPGTYVGPP GGGGPPGTPI MPSPADSTNS SDNIYTMINP VPPGGSRSNF PMGPGSDGPM GGMGMEPHH MNGSLGSGDI DGLPKNSPNN ISGISNPPGT PRDDGELGGN FLHSFQNDNY SPSMTMSV
Source	Baculovirus
Target Names	SSBP3
Protein Names	Recommended name: Single-stranded DNA-binding protein 3 Alternative name(s): Sequence-specific single-stranded-DNA-binding protein
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.