



Recombinant Human Actin-related protein M1 (ARPM1)

Product Code	CSB-YP880135HU
Abbreviation	ARPM1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BYD9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNHCQLPVVI DNGSGMIKAG VAGCREPQFI YPNIIGRAKG QSRAAQGGLE LCVGDQAQDW RSSLFISYPV ERGLITSWED MEIMWKHIYD YNLKPKCDG PVLITEPALN PLANRQQITE MFFEHLGVPA FYMSIQAVLA LFAAGFTTGL VLNSGAGVTQ SVPIFEGYCL PHGVQQDLDA GLDLTNYLMV LMKNHGIMLL SASDRKIVED IKESFCYVAM NYEEEMAKKP DCLEKVYQLP DGKVIQLHDQ LFSCPEALFS PCHMNLEAPG IDKICFSSIM KCDTGLRNSF FSNIILAGGS TSFPGLDKRL VKDIAKVAPA NTAVQVIAPP ERKISVWMGG SILASLSAFQ DMWITAAEFK EVGPNIVHQR CF
Source	Yeast
Target Names	ACTRT3
Protein Names	Recommended name: Actin-related protein M1
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.