



Recombinant Human RanBP-type and C3HC4-type zinc finger-containing protein 1 (RBCK1)

Product Code	CSB-YP880138HU
Abbreviation	RBCK1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BYM8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDEKTKKAAE MALSLTRAVA GGDEQVAMKC AIWLAEQRVP LSVQLKPEVS PTQDIRLWVS VEDAQMHTVT IWLTVRPDMT VASLKDMVFL DYGFPPVLQQ WVIGQRLARD QETLHSHGVR QNGDSAYLYL LSARNTSLNP QELQRERQLR MLEDLGFKDL TLQPRGPLEP GPPKPGVPQE PGRGQPDVAVP EPPVGVWQCP GCTFINKPTR PGCEMCCRAR PEAYQVPASY QPDEEERARL AGEEEALRQY QQRKQQQQEG NYLQHVQLDQ RSLVLNTEPA ECPVCYSVLA PGEAVVLRER LHTFCRECLQ GTIRNSQEA VSCPFIDNTY SCSGKLLERE IKALLTPEDY QRFLDLGSI AENRSAFEYH CKTPDCKGWC FFEDDVNEFT CPVCFHVNCL LCKAIHEQMN CKEYQEDLAL RAQNDVAARQ TTEMLKVMLQ QGEAMRCPQC QIVVQKKDGC DWIRCTVCHT EICWVTKGPR WGPGGPGDTS GGCRCRVNGI PCHPSCQNCH
Source	Yeast
Target Names	RBCK1
Protein Names	Recommended name: RanBP-type and C3HC4-type zinc finger-containing protein 1 EC= 6.3.2.- Alternative name(s): HBV-associated factor 4 Heme-oxidized IRP2 ubiquitin ligase 1 Short name= HOIL-1 Hepatitis B virus X-associated pro
Expression Region	1-510
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is similar to mouse UIP28/UbcM4 interacting protein. Alternative splicing has been observed at this locus, resulting in distinct isoforms.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.