



Recombinant Human Uncharacterized protein KIAA1737 (KIAA1737)

Product Code	CSB-MP880155HU
Abbreviation	KIAA1737
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C0C6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MERKNPSRES PRRLSAKVGK GTEMKKVARQ LGMAAAESDK DSGFSDDGSSE CLSSAEQMES EDMLSALGWS REDRPRQNSK TAKNAFPTLS PMVVMKNVLV KQGSSSSQLQ SWTVQPSFEV ISAQPQLLFL HPPVPSPVSP CHTGEKKSIDS RNYLPILNSY TKIAPHPGKR GLSLGPEEKG TSGVQKKICT ERLGPSLSSS EPTKAGAVPS SPSTPAPPSA KLAEDSALQG VPSLVAGGSP QTLQPVSSSH VAKAPSLTFA SPASPVASD STLHGLESNS PLSPLSANYS SPLWAAEHLC RSPDIFSEQR QSKHRRFQNT LVVLHKSGLL EITLTKELI RQNQATQVEL DQLKEQTQLF IEATKSRAPQ AWAKLQASLT PGSSNTGSDL EAFSDHPAI
Source	Mammalian cell
Target Names	CIPC
Protein Names	Recommended name: Uncharacterized protein KIAA1737
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.