



Recombinant Schizosaccharomyces pombe Zinc-type alcohol dehydrogenase-like protein PB24D3.08c (SPAPB24D3.08c)

Product Code	CSB-BP880165SXV
Abbreviation	SPAPB24D3.08c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C0Y6
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MVSNNAVIK KYLDSTAGYP VIGEHLAFEK REFDLENAQV DEETPVLLKN IYTSVDPYLR MRMQSPKHA S YIPPLELGKP FYNSTVAKVV KSTLDQYKPG MDVVFVSGWE EYTFVSKQAL GFLQPINNPY KLPLIDFVGS LGMPSQTAYC GLKHIGKPKA GETIYISAAS GAVGQMAGQL AKAMGLHVVG SVGSDEKFKI CLDSGYDSVF NYKKESPFKA LPRLCPKGID IYFENVGGET MDAVLENMNL QGRIIFCGAI SQYNNPNPYR VKNLGMVLVK SLTIQGFIVA NILPQYQEY FEEMPKLIAE GKIKYKCDVY DGLESAPEAF IGMLQGKNSG KTIVKIADE
Source	Baculovirus
Target Names	SPAPB24D3.08c
Protein Names	Recommended name: Zinc-type alcohol dehydrogenase-like protein PB24D3.08c EC= 1.-.-.-
Expression Region	1-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.