



Recombinant Schizosaccharomyces pombe Uncharacterized oxidoreductase P32A8.02 (SPAP32A8.02)

Product Code	CSB-MP880173SXV
Abbreviation	SPAP32A8.02
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C1X5
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MESQTESTT VTLTNGMVIP RIGFGAFMLK YNECYGLVTQ ALDSGYRHID TAAVYGNEDI CGKAIVDWCE KNNVKRTDIF LTSKLANCSD YYSTRAAIRS SLHHLGTYID LFLIQSPAGG KKSRIASWKA MEEFVDSGDI RSVGVSNYGV KHLQELYASN PKFYPCVNQI ELHPFLSQDD IVKYCQSHDI AIEAYSPLTH GIRLNDEKLV PIAKKNISV AQLLIRWSLQ KGYIPIIKST KKEHMLSDLD VFNFTIPDDV VQELSSFDEH WHAGTTYDPT VCD
Source	Mammalian cell
Target Names	SPAP32A8.02
Protein Names	Recommended name: Uncharacterized oxidoreductase P32A8.02 EC= 1.-.-.-
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.