



Recombinant Arabidopsis thaliana Glutathione S-transferase F14 (GSTF14)

Product Code	CSB-MP880194DOA
Abbreviation	GSTF14
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C6C8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MADSKMKLHC GFIWGN SAAL FCINEKGLDF ELVFVDWLAG EAKTKTFLST LNPFGEVPVL EDGDLKLFEP KAITRYLAEQ YKDVGTNLLP DDPKKRAIMS MWMEVDSNQF LPIASTLIKE LIINPYQGLA TDDTAVQENK EKLSEVLNIY ETRLGESPYL AGESFSLADL HHLAPIDYLL NTDEEELKNL IYSRPNVA AW VEKMKMRPAW LKTVMKNHI VDLMKQRRLP IKLDSSCHES TVVAQKNAIA IENK
Source	Mammalian cell
Target Names	GSTF14
Protein Names	Recommended name: Glutathione S-transferase F14 Short name= AtGSTF14 EC= 2.5.1.18 Alternative name(s): GST class-phi member 14
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.