



# Recombinant Arabidopsis thaliana AP2/ERF and B3 domain-containing transcription repressor TEM1 (TEM1)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP880196DOA   |
| <b>Abbreviation</b>      | TEM1  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9C6M5  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Arabidopsis thaliana (Mouse-ear cress)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MEYSCVDDSS TTSELSIST TPKPTTTEK KLSSPPATSM RLYRMGSGGS<br>SVVLDSNGV ETESRKLPSS KYKGVVPQPN GRWGAQIYEK HQRVWLGFN<br>EEEEAASSYD IAVRRFRGRD AVTNFKSQVD GNDAESAFLD AHSKAEIVDM<br>LRKHTYADEF EQSRRKFVNG DGKRSGLETA TYGNDAVLRA REVLFEKTVT<br>PSDVGKLNRL VIPKQHAEKH FPLPAMTTAM GMNPSPTKGV LINLEDRTGK<br>VWRFVRSYWN SSQSYVLTGK WSRFVKEKNL RAGDVVCFER STGPDRQLYI<br>HWKVRSSPVQ TVVRLFGVNI FNVSNEKPND VAVECVGKKR SREDDLFLSLG<br>CSKKQAIINI L |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | TEM1  |
| <b>Protein Names</b>     | Recommended name: AP2/ERF and B3 domain-containing transcription repressor TEM1<br>Alternative name(s): Protein TEMPRANILLO 1 RAV1-like ethylene-responsive transcription factor TEM1   |
| <b>Expression Region</b> | 1-361   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                     |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.