



Recombinant Arabidopsis thaliana Agamous-like MADS-box protein AGL92 (AGL92)

Product Code	CSB-BP880198DOA
Abbreviation	AGL92
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C6V4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MRTKTKLVLI PDRHFRRATF RKRNAGIRKK LHELTTLCDI KACAVIYSPF ENPTVWPSTE GVQEVISEFM EKPATERSKT MMSHETFLRD QITKEQNKLE SLRRENRETQ LKHFMFDCVG GKMSEQQYGA RDLQDLSLFT DQYLNQLNAR KKFLTEYGES SSSVPPLFDV AGANPPVVAD QAAVTVPPLF AVAGANLPVV ADQAAVTVPP LFAVAGANLP VVADQAAVNV PTGFHNMNVN QNQYEPVQPY VPTGFSDHIQ YQNMNFNQNG QEPVHYQALA VAGAGLPMTQ NQYEPVHYQS LAVAGGGLPM SQLQYEPVQP YIPTVFSDNV QYQHMNLYQN QQEPVHYQAL GVAGAGLPMN QNQYEPVQPY VPTGFSDHFQ FENMNLNQNQ QEPVQYQAPV DFNHQIQQGN YDMNLNQNMK HAHIPFMDGN YYNYHQPTV GLTSTGHMPS TTTTTTNNNN NNNV
Source	Baculovirus
Target Names	AGL92
Protein Names	Recommended name: Agamous-like MADS-box protein AGL92
Expression Region	1-464
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.