



Recombinant Arabidopsis thaliana 60S ribosomal protein L6-2 (RPL6B)

Product Code	CSB-MP880219DOA
Abbreviation	RPL6B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C9C6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MPAKQRTAKV NRNPDLIRGV GKYSRSQMYH KRGLWAIKAK NGGVFPRHDA KSKVDAPVEK PPKFYPAEDV KKPLPNRRTA KPAKLRASIT PGTVLILAG RFK GKRVVFL KQLASGLLLV TGPFKINGVP LRRVNQAYVI GTSTKVDISG VTLDKFDDKY FGKVAEKKKK KTEGEFFAE KEEKKEIPQG KKDDQKAVDA ALIKAIEAVP ELKTYLGARF SLKQGMKPHE LVF
Source	Mammalian cell
Target Names	RPL6B
Protein Names	Recommended name: 60S ribosomal protein L6-2
Expression Region	1-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.