



Recombinant Arabidopsis thaliana TRM112-like protein At1g78190 (At1g78190)

Product Code	CSB-BP880223DOA
Abbreviation	At1g78190
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C9R3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MRLIVHNMLS CNIKGVVNKF PLRIEAEKVT VKEVDFNPDF LRYMFAKIDW KALVDGARSM EYTELPDNP DTTTLESDT FLRKFHALL ELHLEEGSLV CPETGRKFSV SKGIPNMLLH EDEV
Source	Baculovirus
Target Names	TRM112B
Protein Names	Recommended name: TRM112-like protein At1g78190
Expression Region	1-124
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.