



Recombinant Arabidopsis thaliana Serine carboxypeptidase-like 5 (SCPL5)

Product Code	CSB-BP880236DOA
Abbreviation	SCPL5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9CAU2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	AS IVKFLPGFEG SLPFELETGY IGIGEEEEVQ LFYYFIKSER NPKEDPLLLW LSGGPGCSSI SGLLFENGPL AMKLDVYNGT LPSLVPTTYS WTKTSSMIFL DQPVGTFGSY SRTQQYNKPS DSGEAKRIHE FLQKWLKSHQ EFSSNPFYVA GDSYSGMVVP ATVQEISKGN YQCCSPPINL QGYVLGNPIT EHAIDYNYRI PFAHGMALIS DELYESLKRV CKGEYVDPDRD TECLKLVEEF SKCTKGVQCQE VVIKPLCVTE TPNCYIYRYL LTTYWVNDVN VRKALQINKE SIGEWVRCYF GIPYTHDIKS SVPYHMNSI NGYRSLIYSG DHDLNVPFLA TQAWVRS LNY SIIDNWRPWW IKDQIGGYTK TYANKMTFAT VRGGGHTAEY KPYETYIMFH RWINGQPL
Source	Baculovirus
Target Names	SCPL5
Protein Names	Recommended name: Serine carboxypeptidase-like 5 EC= 3.4.16.-
Expression Region	29-438
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.