



# Recombinant Mouse Keratin, type I cytoskeletal 20 (Krt20)

<b>Product Code</b>	CSB-MP880486MO
<b>Abbreviation</b>	Krt20
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9D312
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDFSRQSFHR SLSSSSQGPA LMSGSLYRK GTVQRLGAAP SVYGGAGGHG TRISVSKAVM SYGGDLSNGS DLFGGNGKLA MQNLNDRLAN YLEKVRSEIQ SNSRLEAIQK QWYETNAPST IRDYSSYYAQ IKELQNQVKD AQVQNAQCVL RIDNAKLAEE DFRLKFETER GMRIAVEADL QGLSKVYDNL TLQKTDLEIQ IEELNKDLAL LKKEHQEEVE VLRRLGNNV NVEVDAAPGL NLGEIMNEMR QRYEVLAQKN LQEAKEQFER QSQTLQQQVT VNTEELKGFV VQVTELRRTY QNLEIELQSH LSMKESLERN LEDVKARYAS QLAAIQEMLS SLEAQLMQIR SDTERQNQEH NILLDIKTRL EQEIATYRRL LEGEDIKTTTE YQLSTLEMKD IKKTRKIKTV VEEVVDGKVV SSEVKEIEES V
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Krt20
<b>Protein Names</b>	Recommended name: Keratin, type I cytoskeletal 20 Alternative name(s): Cytokeratin-20 Short name= CK-20 Keratin-20 Short name= K20
<b>Expression Region</b>	1-431
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein is a member of the keratin family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. The type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. This cytokeratin is a major cellular protein of mature enterocytes and goblet cells and is specifically expressed in the gastric and intestinal mucosa. The type I cytokeratin genes are clustered in a region of chromosome 17q12-q21.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.