



Recombinant Mouse Sterol 26-hydroxylase, mitochondrial (Cyp27a1)

| | |
|--------------------------|---|
| Product Code | CSB-YP880560MO |
| Abbreviation | Cyp27a1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9DBG1 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | TIPAAALQA QESTEGPGTG QDRPRLRSPA ELPGTGTLQF LFQLFLQGYV LHLPLDLQVLN KTKYGPMTWT SFGTYTNVNL ASAPLLEQVM RQEGKYPIRD HMDQWKDHRD HKGLTYGIFI AQGEQWYHLR QALKQRLLKP DEAAALYTDAL NEVISDFITR LDQVRAESES GDQVPDMAHL LYHLALEAIT YILFEKRIGC LKPSIPEDTA AFIRSVAIMF QNSVYITFLP KWTRPLLPFW KRYLNGWDNI FSFGKKLIDE KVQELKAQLQ ETGPDGVRVS GYLHFLTNE LLSTQETIGT FPELLLAGVD TTSNTLTWAL YHLSKSPEIQ EALHKEVTGV VPFGKVPQHK DFAHMPLLKA VIKETRLRYP VVPTNSRIIT EKETEINGFL FPKNTQFVLC HYVVS RDPSV FPEPNSFQPH RWLRKKEADN PGILHPFGSV PFGYGVR SCL GRRIAELEMQ LMLSRLVQKY EIALAPGMGE VKTVSRIVLV PSKKVRLHFL QRQ |
| Source | Yeast |
| Target Names | Cyp27a1 |
| Protein Names | Recommended name: Sterol 26-hydroxylase, mitochondrial EC= 1.14.13.15 Alternative name(s): 5-beta-cholestane-3-alpha,7-alpha,12-alpha-triol 27-hydroxylase Cytochrome P-450C27/25 Cytochrome P450 27 Sterol 27-hydroxylase Vitam |
| Expression Region | 33-533 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Target Details | This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This mitochondrial protein oxidizes cholesterol intermediates as part of the bile synthesis pathway. Since the conversion of cholesterol to bile acids is |



the major route for removing cholesterol from the body, this protein is important for overall cholesterol homeostasis. Mutations in this gene cause cerebrotendinous xanthomatosis, a rare autosomal recessive lipid storage disease.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.