



Recombinant Rat Platelet-derived growth factor D (Pdgfd)

Product Code	CSB-MP880633RA
Abbreviation	Pdgfd
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9EQT1
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	TPQSASI KALRNANLRR DESNHLTDLY RRDENIRVTG TGHVQSPRFP NSYPRNLLLT WRLHSQEKTR IQLAFDHQFG LEEAENDICR YDFVEVEDVS ESSTVVRGRW CGHKEIPPRI TSRTNQIKIT FQSDDYFVAK PGFKIYYSFV EDFQPEAASE INWESVTSSF SGVSYHSPSV MDSTLTADAL DKAIAEFDTV EDLLKYFNPA SWQDDLENLY MDTPRYRGRS YHERKSKVDL DRLNDDVKRY SCTPRNHSVN LREELKLTNA VFFPRCLLVQ RCGGNCGCGT LNWKSCTCSS GKTVKKYHEV LKFEPGHFKR RGKAKNMALV DIQLDHHERC DCICSSRPPR
Source	Mammalian cell
Target Names	Pdgfd
Protein Names	Recommended name: Platelet-derived growth factor D Short name= PDGF-D Alternative name(s): Iris-expressed growth factor Spinal cord-derived growth factor B Short name= SCDGF-B Cleaved into the following 2 chains: 1. Platelet-
Expression Region	24-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a member of the platelet-derived growth factor family. The four members of this family are mitogenic factors for cells of mesenchymal origin and are characterized by a core motif of eight cysteines, seven of which are found in this factor. This gene product only forms homodimers and, therefore, does not dimerize with the other three family members. It differs from alpha and beta members of this family in having an unusual N-terminal domain, the CUB domain. Two splice variants have been identified for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.