



Recombinant Mouse Tumor necrosis factor receptor superfamily member 23 (Tnfrsf23)

Product Code	CSB-BP880638MO
Abbreviation	Tnfrsf23
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9ER63
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	M PESYSFNCPD GEYQSNDVCC KTCPSGTFVK APCKIPHTQG QCEKCHPGTF TGKDNLHDC ELCSTCDKDQ NMVADCSATS DRKCECQIGL YYYPKFPES CRPCTKCPQG IPVLQECNST ANTVK
Source	Baculovirus
Target Names	Tnfrsf23
Protein Names	Recommended name: Tumor necrosis factor receptor superfamily member 23 Alternative name(s): Decoy TRAIL receptor 1 TNF receptor family member SOB TNF receptor homolog 1 Tumor necrosis factor receptor p60 homolog 1
Expression Region	30-155
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.