



Recombinant Arabidopsis thaliana Putative F-box protein At5g38390 (At5g38390)

Product Code	CSB-MP880713DOA
Abbreviation	At5g38390
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FF23
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDLLSNLPYE LLCHILSFLT TKEAALTSVL SKRWRNLI AF VPNVDIVDHD IRELFMDFVD RVLALQGN SP INKFSLDCSG VSDRVDCWI QNVMVRGVSE LNL SIFVDSV FEDDYNLFPK VFENKKLVKL GLSYISWLDG SIFLPMLKTL VLESVLLSVE KFEILLHALP ALEELVMNNI YWKERDVVDL IVSSESLKTL TINFIVCTHT LSFDTPI LAY LSYSGYVVD D YAEAKMENLF EARISLLVVE GDISRVRALI NNDLLEDDEY DVLQFENVWK LMNGIRNIRC LYLSPNTLEV LSLCCESMPV FKNLKSLSIK SAENRGWQAM PVLLRNCPHL ETLVLEGLLH HVTDKCGDAC DCVSREDKGR SLTSCPVKVL EIKGFQGTK EMNMIKHFL E YFPCLKELKM YMEENDPAQL RVPEVITGEM MERYKKSSSC NVQLLVGGKL YKKWTP
Source	Mammalian cell
Target Names	At5g38390
Protein Names	Recommended name: Putative F-box protein At5g38390
Expression Region	1-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.