



Recombinant Arabidopsis thaliana Ethylene-responsive transcription factor ABR1 (ABR1)

Product Code	CSB-YP880727DOA
Abbreviation	ABR1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FGF8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MCVLKVANQE DNVGKKAESI RDDDHRTLSE IDQWLYLFAA EDDHHRHSFP TQQPPSSSS SSLISGFSRE MEMSAIVSAL THVVAGNVPQ HQQGGGEGSG EGTSNSSSSS GQKRREVEE GGAKAVKAAN TLTVDQYFSG GSSTSKVREA SSNMSGPGPT YEYTTTATAS SETSSFSGDQ PRRRYRGVRQ RPWGKWAAEI RDPFKAARVW LGTFDNAESA ARAYDEAALR FRGNKAKLNF PENVKLVRPA STEAQPVHQT AAQRPTQSRN SGSTTTLLPI RPASNQSVHS QPLMQSYNLS YSEMARQQQQ FQQHHQQSLD LYDQMSFPLR FGHTGGSMQ STSSSSSHSR PLFSPAAVQP PPESASETGY LQDIQWPSDK TSNNYNNSPS S
Source	Yeast
Target Names	ABR1
Protein Names	Recommended name: Ethylene-responsive transcription factor ABR1 Alternative name(s): Protein ABA REPRESSOR 1
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.