



Recombinant Arabidopsis thaliana Putative F-box/FBD/LRR-repeat protein At5g52460 (EDA41)

| | |
|--------------------------|--|
| Product Code | CSB-BP880738DOA |
| Abbreviation | EDA41 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9FHB5 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MAERCQRIDE EVVASRDEIS SLPDDLLIQI LLLVPIKDAV GTMILSKRWR YVWTLLEPKLE YSDPGDECES VWKFLEKLE LHKAPFLETL CIQLDVGKLV ANAVDRFVRK LELELHWTAE PTSLPKSLYT CKTLVELTSLS DKIIVDVPSS VCLPSLNILR LFYVVFKDEN SLERLISSCS VLARSKVGLA ISDNRDYRFA KVHKLRLKLP HIDVVCHTDD KFLKAISLVT YLVLPLEDPM ALDFRGFTFT RLGKLVICPH GHLWLDIPL ILNNSPKLRF LAITSVMDID PEISHCHGTN QGTVPRCLSA HLDEEFVWHG YRGNEEETQL IRYIFANAKC LKKREISTFH LEEREVIETV LKSMPRVSTT STLVFK |
| Source | Baculovirus |
| Target Names | EDA41 |
| Protein Names | Recommended name: Putative F-box/FBD/LRR-repeat protein At5g52460 Alternative name(s): Protein EMBRYO SAC DEVELOPMENT ARREST 41 |
| Expression Region | 1-376 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |