



Recombinant Arabidopsis thaliana Putative F-box/FBD/LRR-repeat protein At5g52460 (EDA41)

Product Code	CSB-EP880738DOA
Abbreviation	EDA41
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FHB5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MAERCQRIDE EVVASRDEIS SLPDDLLIQI LLLVPIKDAV GTMILSKRWR YVWTLTPKLE YSDPGDECES VWKFLEKLE LHKAPFLETL CIQLDVGKLV ANAVDRFVRK LELELHWTAE PTSLPKSLYT CKTLVELTLS DKIIVDVPSS VCLPSLNILR LFYVVFVDEN SLERLISSCS VLARSKVGLA ISDNRDYRFA KVHKLRLKLP HIDVVCHTDD KFLKAISLVT YLVLPLEDPM ALDFRGFTFT RLGKLVICPH GHLWLDIPL ILNNSPKLRF LAITSVMDID PEISHCHGTN QGTVPRLSA HLDEEFVWHG YRGNEEETQL IRYIFANAKC LKKREISTFH LEEREVIETV LKSMPRVSTT STLVEK
Source	E.coli
Target Names	EDA41
Protein Names	Recommended name: Putative F-box/FBD/LRR-repeat protein At5g52460 Alternative name(s): Protein EMBRYO SAC DEVELOPMENT ARREST 41
Expression Region	1-376
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.