



Recombinant Arabidopsis thaliana Ubiquitin-conjugating enzyme E2 27 (UBC27)

Product Code	CSB-BP880746DOA
Abbreviation	UBC27
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FI61
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MIDFSRIQKE LQDCERNQDS SGIRVCPKSD NLTRLTG TIP GPIGTPYEGG TFQIDITMPD GYPFEPKMQ FSTKVWHPNI SSQSGAICLD ILKDQWSPAL TLKTALVSIQ ALLSAPEPKD PQDAVVAEQY MKNYQVFVST ARYWTETFAK KSSLEEKVKR LVEMGFGDAQ VRSAIESSGG DENLALEKLC SA
Source	Baculovirus
Target Names	UBC27
Protein Names	Recommended name: Ubiquitin-conjugating enzyme E2 27 EC= 6.3.2.19 Alternative name(s): Ubiquitin carrier protein 27
Expression Region	1-192
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.