



Recombinant Arabidopsis thaliana UDP-glycosyltransferase 76C3 (UGT76C3)

Product Code	CSB-BP880747DOA
Abbreviation	UGT76C3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FI96
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDKSNGLRVI LFPLPLQGCI NPMIQLAKIL HSRGFSITVI HTRFNAPKAS NHPLFTFLQI PDGLSETETR THDITLLLTL LNRSCESPFR ECLKLLQSA DSETGEEKQR ISCLIDDSGW IFTQPVAQSF NLPRLVLNTY KVSFFRDHFV LPQLRREMYL PLQDSEQGDD PVEEFPPLRK KDLLQILDQE SEQLDSYSNM ILETTKASSG LIFVSTCEEL DQDSLSQARE DYQVPIFTIG PSHSYFPGSS SSLFTVDETC IPWLDKQEDK SVIYVSFGSI STIGAEFME IAWALRNSDQ PFLWVVRGGS VVHGAEWIEQ LHEKGKIVNW APQQEVLKHQ AIGGFLTHNG WNSTVESVFE GVPMICMPFV WDQLLNARFV SDVWMVGLHL EGRIERNVIE GMIRRLFSET EGKAIRERME ILKENVGRSV KPKGSAYRSL QHLIDYITYF
Source	Baculovirus
Target Names	UGT76C3
Protein Names	Recommended name: UDP-glycosyltransferase 76C3 EC= 2.4.1.-
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.