



Recombinant Arabidopsis thaliana Germin-like protein subfamily 1 member 17 (At5g39150)

Product Code	CSB-MP880748DOA
Abbreviation	At5g39150
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FIC6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	YDPSPLQDF CVAIDDPKNG VFVNGKFCKD PKQAKAEDFF SSGLNQAGIT NNKVQSNVTT VNVDQIPGLN TLGISLVRID YAPYGQNPPH THPRATEILV LVEGTLYVGF VSSNQDNNRL FAKVLNPGDV FVFPIGMIHF QVNIGKTPAV AFAGLSSQNA GVITIADTVF GSTPPINPDI LAQAFQLDVN VVKDLEAKFK N
Source	Mammalian cell
Target Names	At5g39150
Protein Names	Recommended name: Germin-like protein subfamily 1 member 17
Expression Region	22-221
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.