



Recombinant Arabidopsis thaliana F-box protein At5g18160 (At5g18160)

Product Code	CSB-BP880769DOA
Abbreviation	At5g18160
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FK54
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDKQDEKKQG TTKSSSTLTT RCSHGNIHSQ SNSIPLDITI EILSRPAAKS IVRSRSVSKL WSSITTTPEF IKHRSKKTSP PCVLLIFRKH DKLIVFSSPQ HQNTYSHVQD YHIEIPKNGF IRRLDVHGL ICLEGSKQLV ICNPTLKRFF PLPEPQGTGD EYNVGGFLGY EPIEGKYKAL CIVRGWNTQV LTLEIQESWR VTKPGYTHWP TKDTGRCING VIYYKAIIFD RVPRHVIFGF DLRYEETFHI EFPMRNYDRF LMVSYEGRLA LISSTSSVVE IWSLEDAGNR KWSYEQFHLC LPPNTSLKGV IDAGELIYTG FSLNRSFCVV YFDPKKNIR ETKFQGIADN QIWQPDLRGF DLVNDFYVLP NHIESFISF
Source	Baculovirus
Target Names	At5g18160
Protein Names	Recommended name: F-box protein At5g18160
Expression Region	1-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.