



Recombinant Arabidopsis thaliana Serine/threonine-protein phosphatase 7 (PP7)

Product Code	CSB-BP880803DOA
Abbreviation	PP7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FN02
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	METVPPSPIT WPDGGALTND WWHGLMSCFE WSSWNLPPSQ LPSLLPVNVF DSLVLTAHKI LHKERNCVHI DDLDSVSNVV VVGDIHGQLH DLLFLLKDTG FPCQNR CYVF NGDYVDRGAW GLETFLVLLS WKVLM PDRVY LLRGNHESKY CTSMYGF EKE VLTKYGD K GK HVYRKCLGCF EGLPLASIIS GRVYTAHGGL FRSPVLPKRT TRGKKNRRVV LLEPEPSSMK LGTLDELMQA RRSVLDPPWE GSNLIPGDVL WSDPSMTPGL SPNEQRGIGL LWGPDC TEDF LKKYELKLI RSHEGPDARE KRTGLGGMDN GYTIDHNVES GKLITIFSAP DYPQFQATEE RYKNKGAYII LQAPDFSDPQ FHSFEAVKPR PKAHPYYDFE NVIDSDDEMD KSAMDTNNEQ PNS
Source	Baculovirus
Target Names	PP7
Protein Names	Recommended name: Serine/threonine-protein phosphatase 7 EC= 3.1.3.16
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.