



Recombinant Arabidopsis thaliana UDP-glycosyltransferase 79B9 (UGT79B9)

Product Code	CSB-EP880804DOA
Abbreviation	UGT79B9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FN28
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGQNFHAFMF PWFAFGHMTP YLHLANKLAA KGHRVTFLLP KKAQKQLEHH NLFPDRIIFH SLTIPHVDGL PAGAETASDI PISLGKFLTA AMDLTRDQVE AAVRALRPDL IFFDTAYWVP EMAKEHRVKS VIYFVISANS IAHELVPGGE LGVPPPGYPS SKVLYRGHDA HALLTFSIFY ERLHYRITTG LKNCDFISIR TCKEIEGKFC DYIERQYQRK VLLTGPMLE PDNSRPLEDR WNHWLNQFKP GSVIYCALGS QITLEKDQFQ ELCLGMELTG LPFLVAVKPP KGAKTIQEAL PEGFEERVKN HGVVWGEWVQ QPLILAHPSV GCFVTHCGFG SMWESLVSDC QIVLLPYLCD QILNTRLMSE ELEVSVEVKR EETGWFSKES LSVAITSVMD KDSELGNLVR RNHAKLKEVL VSPGLLTGYT DEFVETLQNI VNDTNLE
Source	E.coli
Target Names	UGT79B9
Protein Names	Recommended name: UDP-glycosyltransferase 79B9 EC= 2.4.1.-
Expression Region	1-447
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.